

Invacare® Top End® Terminator™ Jr. Everyday Wheelchair

Price List and Order Form Effective July 1, 2003 Revised October 25, 2004



Yes, you can.™

For accuracy of ordering, please fax to Top End at 727-522-1007.

For questions call 1800-532-8677 or 727-522-8677.

To order Literature: 800-828-6282

www.invacare.com

Provider Name _____ Date _____ Top End Job Number _____
 Account # _____ P.O. # _____ Phone # _____ Fax # _____
 Address _____ City _____ State _____ Zip Code _____
 Customer Name _____ Phone No _____
 Address _____ City _____ State _____ Zip Code _____
Ship order to: **Provider** **Customer** **Order will ship to Provider unless otherwise noted.**

Chair Type
 TERJR Invacare Top End Terminator Jr. Everyday\$2,137.00 _____
 Weight limit - 150 lbs

Seat width: 10" 11" 12" 13" 14" 15" 16"

Adjustables seat depth 10" - 18" :
 10" 11" 12" 13" 14" 15" 16" 17" 18"

Adjustable frame height choose high or low:
 (See seat-to-floor matrix on back)
 High seat-to-floor 16" - 20"
 Low seat-to-floor 13" - 18"
Set at: front _____ **rear** _____

Seat: Rigid seat pan w/adjustable angle 0 - 2"
 TERJRI Cushion 2" soft.....75.00 _____

Back: Sportster (w/o push handles 10" - 16")standard
 10° Bend back w/push handles (17" - 25" tall).....**No charge**

Back Accessories:
 TERJR2 Adjustable height removable push handles (17" - 21")200.00 _____
 TERJR3 Stroller handles (height from seat: 24" - 31")100.00 _____

Camber:
 Adjustable A4 camber system nut & camber insert release
²Camber angle: widths 12"-16"
 (choose two) 0 3 6 9 12
 Camber angle: widths 10" & 11"
 (choose one) 0 3 6 9 12
 TERJR4 Quick release lever.....30.00 _____

Wheel locks High (push to lock):
 TERJR5 Wheel lock extension (8")30.00 _____

Back upholstery: Adjustable tension w/center and 2 flaps
 TERJR7 Custom manual adjustable tension upholstery.....100.00 _____
 TERJR6 Back pack50.00 _____

Rear wheel size:
⁴ 18" 20" 22" 24"

Wheel choice:
 Spoke.....standard
 Composite.....**No charge**

Tires - 13/8" pneumatic treaded:
 TERJR8 Pneumatic w/FF inserts.....125.00

Pushrims - anodized:
 TERJR9 Coated pushrims100.00

Spoke guards:
 TERJR10 Spoke protectors (clear).....75.00 _____
 TERJRI1 Smiley face
 20" 22" 24"150.00 _____
 Flower (22" only)150.00 _____
 Spiderweb (22" only)150.00 _____

Front wheels size semi-pneumatic:
 4" 5" 6"⁵
 TERJR12 4" Lighted casters70.00 _____

Footboard flip-up height/angle adjustable:
⁷ Flare front caster housing.....**No charge**

Arm options:
 TERJR20 Cantilever arms150.00 _____
⁸ **TERJR21** T-arms adjustable height desk ...175.00 _____
⁹ **TERJR22** Removable side guard/fender ..175.00 _____
 TERJR14 Removable side guards125.00 _____

Safety options:
 Rear anti-tippersstandard
 TERJR15 Crutch holder w/straps.....125.00 _____
 TERJR18 Hook & loop seat positioning strap..40.00 _____
 TERJR19 Auto style buckle seat positioning strap50.00 _____

Special Instructions

Frame color

CONTEMPORARY COLORS

Wet Black Black Prism Midnight Blue
 Silver Metallic Silver Vein

WAY COOL COLORS

Electric Blue Lolly Pop Blue Cosmic Blue
 Celestial Purple Vivid Violet Electric Teal
 Emerald Green Shamrock Green Grasshopper Green

TOO HOT COLORS

Tangerine Red Electric Red Grape Madness
 Sunny Yellow Sunkissed Orange
 Fruit Fusion Bubble Gum Pink

- FOOTNOTES**
1. Must be ordered with 10° back with push handles.
 2. Only one camber available with narrow widths of 10" and 11".
 3. Adjustable tension by straps.
 4. 18" not available in spoke.
 5. 6" not available w/18" or 20" rear wheel.
 6. Footboard range from seat:
High frame 2.5" - 12.75". Low frame 2.5" - 8.5".
 7. Recommended only for footboard settings less than 5" from floor.
 8. Recommended for use with lap tray.
 9. Not available with T-arms.

Invacare® Top End® Terminator™ Jr. Everyday Wheelchair

Measuring Guidelines

Low Frame Matrix			
Rear Wheel Wheel Size	Mounting Position	Front Seat Height	Rear Seat Height
18"	Upper	13"	11"
	Center	14"	12"
	Lower	15"	13"
20"	Upper	14"	12"
	Center	15"	13"
	Lower	16"	14"
22"	Upper	15"	13"
	Center	16"	14"
	Lower	17"	15"
24"	Upper	16"	14"
	Center	17"	15"
	Lower	18"	16"

High Frame Matrix			
Rear Wheel Wheel Size	Mounting Position	Front Seat Height	Rear Seat Height
20"	Upper	16"	14"
	Center	17"	15"
	Lower	18"	16"
22"	Upper	17"	15"
	Center	18"	16"
	Lower	19"	17"
24"	Upper	18"	16"
	Center	19"	17"
	Lower	20"	18"

Back
<ul style="list-style-type: none"> ♦ Back rest tubing is 1" in diameter for seating systems.

Push Handles/Stroller Handles Option
<ul style="list-style-type: none"> ♦ 10° Bend back is equipped with push handles that extend out of the top of the back canes. They adjust in height from 17" - 25" from seat pan. These push handles can be removed if needed to improve propulsion. 1" caps are available upon request. ♦ Stroller handles if ordered must be ordered in combination with the 10° back with push handles. Height from seat is 24" - 31" ♦ Adjustable height push handles bolt to the backrest spreader bar and are recommended for customers who do not use push handles on a regular basis. The height adjusts from 17" - 21" above seat pan.

Flip Up Footboard
<ul style="list-style-type: none"> ♦ Height adjustment: High frame 2.5" - 12.75" from seat pan Low frame 2.5" - 8.5" from seat pan ♦ A flared front caster housing is recommended for customers expected to grow beyond ranges. Longer footboard attachments are available upon request.

Arm Options
<ul style="list-style-type: none"> ♦ Adjustable height padded T-arms are recommended for use with a lap tray. ♦ Cantilever arms are adjustable 7" - 12" high non-locking tubular construction and foam padded.

Specifications	
♦ Frame style	Rigid
♦ Seat widths	10" - 16"
♦ Seat depths	10" - 18"
♦ Adjustable height back	10" - 16"
♦ Adjustable angle folding back	-4°, 0°, 4°, 8°, 12°, 16°
♦ Camber angles	0°, 3°, 6°, 9°, 12°
♦ Tubing O.D.	1"
♦ Weight limit.	150 lbs.
♦ Weight of chair	from 27 lbs.
♦ Warranty	Lifetime on frame and crossmembers

Minimum Overall Width
Seat width plus 9" Measurement taken from outside of handrim to outside of handrim with 0° camber.

Weight Limitation 150 Pounds
Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.
NOTE: All specifications and dimensions are approximate.

Invacare Corporation www.invacare.com



INVACARE TOP END
 4501 63rd Circle North
 Pinellas Park, FL USA
 33781-5914
 800-532-8677

Invacare, Top End and the medallion design are registered trademarks and Terminator and Yes, You can. are trademarks of Invacare Corporation.
 Specifications and prices are subject to change without notification.
 ©2004 Invacare Corporation

Form No. 00-289 Rev 10-25-04
 P2000_289F